

When the job demands contaminant-free transportation of refrigerants, oxygen and other gases, CerroPure® nitrogen charged tube is a must. Specially cleaned for refrigeration, medical and other specialized equipment applications, CerroPure gives you what no other can - Cerro's reputation for quality, service, trust and dependability.

To ensure exceptional, on-the-job cleanliness, CerroPure is purged with dry nitrogen and the ends sealed with rubber plugs to maintain a positive nitrogen charge. In addition, CerroPure tube is clearly marked "N2 Gas" for positive product identification, meets the ASTM B280 standard for ACR tubing, and, is available in standard 10- and 20-foot lengths, 3/8" to 3-1/8" diameter.

For all your copper tube requirements, contact or ask for Cerro. A world leader in copper tube with a continued commitment to uncompromising quality and new service offerings to the markets we support.

WIN! Go to www.cerroflow.com/clean to learn more about Cerro's nitrogen charged tube and register to win a \$50 gift certificate of your choice.



f. 618-332-0108 www.cerroflow.com







